Rotor Safety

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Español

Italiano

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Every Beckman Coulter rotor is manufactured and tested with your safety in mind; however, improper use can result in safety hazards. **Please read this booklet.** It contains information critical to the safe operation and maintenance of rotors and accessories.

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Introduction

This safety booklet summarizes information basic to the safe operation of the rotor described in the instruction manual shipped with the rotor. The international symbol displayed is a reminder that all safety instructions should be read and understood before installation, operation, maintenance, or repair of the rotor is attempted. When you see the international safety symbol above in rotor and centrifuge manuals, pay special attention to the safety information presented. Observing safety precautions will help to prevent actions that could cause personal injury or rotor damage, or that could adversely affect rotor performance. We recommend that you read this safety booklet, the centrifuge manual, and the applicable rotor manual before operating your centrifuge.

Mechanical Safety

- Use only Beckman Coulter rotors designed for use in the centrifuge you are operating. Each Beckman Coulter rotor is developed, manufactured, and tested for safety and reliability as part of a centrifuge/rotor system. The safety and reliability of a rotor cannot be assured if used in a non-Beckman Coulter centrifuge or in a Beckman Coulter centrifuge that has been modified without Beckman Coulter's approval.
- Use only components and accessories designed for use in the rotor being centrifuged (refer to the applicable rotor manual). *The safety of rotor components and accessories made by other manufacturers cannot be ascertained by Beckman Coulter. Use of other manufacturers' components or accessories (other than those listed in the rotor manual) in Beckman Coulter centrifuge rotors may void the rotor warranty and should be prohibited by your laboratory safety officer.*
- The strength of glass and plastic containers can vary between lots, and will depend on handling and usage. We highly recommend that you pretest labware in the rotor (using buffer or

gradient of equivalent density to the intended sample solution) to determine optimal operating conditions. Scratches (even microscopic ones) significantly weaken glass tubes.

- To prevent plug damage on vertical tube or near-vertical tube rotors, do not use rotor plugs or spacers in empty cavities. Leave unused cavities completely empty.
- Ensure that the buckets (or carriers) are properly hooked to the crossbars or pivot pins of swinging bucket rotors. A mishooked bucket cannot swing freely to a horizontal position during centrifugation and could cause a rotor mishap. Attach all buckets, loaded or empty, to the rotor as described in the rotor manual.
- Before starting the centrifuge, make sure that the rotor tie-down device—in centrifuges requiring them—is securely fastened.
- Never exceed the maximum rated speed of the rotor in use. Rotors are designed for use at the speeds indicated (engraved on the rotor and listed in the rotor manual). Also, speed reductions may be required because of weight considerations—the weight of a tube, cap, adapter, or the density of the solution being centrifuged. Be sure to follow the instructions in the applicable rotor manual.
- NEVER attempt to slow or stop a rotor by hand.
- If glass tubes break, remove the glass very carefully from adapters, buckets, or rotor cavities. If all the glass particles are not contained in the bucket or cavity, be careful when examining or cleaning centrifuge gaskets and chambers, as glass particles may be embedded in their surfaces.
- Do not use sharp tools on the rotor that could cause scratches in the rotor surface. Corrosion begins in scratches and may open fissures in the rotor with continued use.

- To prevent premature failures or hazards by detecting stress corrosion, metal fatigue, wear or damage to anodized coatings, and to instruct laboratory personnel in the proper care of rotors, Beckman Coulter offers the Field Rotor Inspection Program (FRIP). This program involves a visit to your laboratory by a specially trained Beckman Coulter representative, who will inspect all your rotors for corrosion or damage. The representative will recommend repair or replacement of at-risk rotors to prevent potential rotor failures. Contact your local Beckman Coulter office to request this service.

Safety Against Risk of Fire

The rotors and accessories are not designed for use with materials capable of developing flammable or explosive vapors. Do not centrifuge such materials in nor handle or store them near the centrifuge.

Chemical and Biological Safety

Normal operation may involve the use of solutions and test samples that are pathogenic, toxic, or radioactive. Do not use such materials in these rotors, however, unless you take *all necessary safety precautions*.

- Observe all cautionary information printed on the original solution containers.

- Handle body fluids with care because they can transmit disease. No known test offers complete assurance that such fluids are free of micro-organisms. Some of the most virulent—Hepatitis (B and C) viruses, HIV (I–V), atypical mycobacteria, and certain systemic fungi—further emphasize the need for aerosol protection. Handle other infectious samples according to good laboratory procedures and methods to prevent spread of disease. Because spills may generate aerosols, observe proper safety precautions for aerosol containment. Do not run toxic, pathogenic, or radioactive materials in this rotor without taking appropriate safety precautions. Biosafe containment should be used when Risk Group II materials (as identified in the World Health Organization *Laboratory Biosafety Manual*) are handled; materials of a higher group require more than one level of protection.
- Dispose of all waste solutions according to appropriate environmental health and safety guidelines.
- If disassembly reveals evidence of leakage, you should assume that some fluid escaped the rotor. Apply all appropriate decontamination procedures to the centrifuge and accessories.

 **NOTE**

Continuous flow, elutriation, and zonal rotors may require additional safety precautions. Carefully observe all cautions and warnings in the applicable rotor manual when using one of these rotors.

It is your responsibility to decontaminate rotors and accessories before requesting service by a Beckman Coulter Field Service representative.

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All rights reserved
Beckman Coulter, Inc. • 250 S. Kraemer Blvd. • Brea, California 92821
©2009 Beckman Coulter, Inc.